

Time Table For Portable Oxygen Cylinders (“E” size)

FLOW:	10 Lpm	6 Lpm	5 Lpm	4 Lpm	3 Lpm	2 Lpm	1 Lpm
Tank pressure	Hr:min						
2200	1:05	1:50	2:15	2:55	3:55	5:55	12:05
2100	1:00	1:45	2:10	2:45	3:45	5:40	11:30
2000	0:55	1:40	2:05	2:35	3:30	5:25	11:00
1900	0:50	1:35	1:55	2:30	3:20	5:05	10:25
1800	0:50	1:30	1:50	2:20	3:10	4:50	9:50
1700	0:45	1:25	1:45	2:10	3:00	4:35	9:15
1600	0:40	1:20	1:35	2:05	2:50	4:15	8:45
1500	0:40	1:15	1:30	1:55	2:35	4:00	8:10
1400	0:35	1:05	1:25	1:45	2:25	3:45	7:35
1300	0:30	1:00	1:15	1:40	2:15	3:25	7:05
1200	0:30	0:55	1:10	1:30	2:05	3:10	6:30
1100	0:25	0:50	1:05	1:20	1:50	2:55	6:00
1000	0:20	0:45	0:55	1:15	1:40	2:35	5:25
900	0:20	0:40	0:50	1:05	1:30	2:20	4:50
800	0:15	0:35	0:45	0:55	1:20	2:05	4:15
700	0:10	0:30	0:35	0:50	1:10	1:45	3:45
600	0:10	0:25	0:30	0:40	0:55	1:30	3:10
500	0:05	0:20	0:25	0:30	0:45	1:15	2:35
400	0:05	0:10	0:15	0:25	0:35	0:55	2:05
300	0:03	0:05	0:10	0:15	0:25	0:40	1:30

Refill Tank when in Gray (Red) Zone 500psi